

# BUCKET LINERS



## BUCKET LINERS

### DESCRIPTION

Pharmaclean<sup>®</sup> LDPE bucket liners for cleanroom, thickness 100 µm, dimensions 800x500 mm. Heat-sealing closure. Opening on the long side (800 mm). Produced and packed in cleanroom ISO 7 according to GMP regulations. They can be used to create a barrier between the sanitizer and the bucket. The bucket liner implementation reduces the possible cross contamination of the disinfectant solution. They are suitable for 15 L and 22 L capacity.

The product is available also in sterile version (Gamma irradiation). Each sterile pack has red dot indicator, code, lot number and expiry date. For each supply of sterile product we provide related irradiation certificate.

# BUCKET LINERS

GENERAL PROPERTIES	
film material	LDPE
film thickness	100 µm
Colour	Transparent
Printing	None
Welding	Bottom side
Additives	Slip, antiblock
Compliance	Foodgrade Declaration EU 10/2001
Storage conditions	The bags must be protected from light, preferably in their original packaging at a storage temperature of ideally 15°C - 30°C and relative humidity of 40-65% RH, stored and protected against moisture and exposure to dirt and other foreign matters

PRODUCTS		
CODE	DESCRIPTION	PACKAGING
PH10118	Pharmaclean LDPE bucket liner, dim. 800 x 500 mm, thickness 100 µm	Pack 25 pieces in double bag
PH10118S	Pharmaclean LDPE sterile bucket liner, dim. 800 x 500 mm, thickness 100 µm	Pack 25 pieces in double bag

COMPANY WITH  
QUALITY SYSTEM  
CERTIFIED BY DNV GL  
= ISO 9001:2015 =