



# PHARMACLEAN® CLEANROOM MATS

#### **DESCRIPTION**

Pharmaclean® tear-off decontamination mats with antimicrobial additive for entry into contamination-controlled areas. The additive is contained in the adhesive part of the mat, consisting of a low-density polyethylene film, and remains on it for long time, ensuring a substantial reduction in cross-contamination. Numbered corners to prevent tearing multiple layers at once. Available in blue, white and grey color and in packs of 8 mats of 30 strips each, or 4 mats of 60 strips each.





Pharmaclean® tear-off mats are the first contamination control product containing renewable plant materials, offering an environmentally friendly alternative to support sustainability.

#### A real commitment

Pharmaclean® mats are made of EnviroTack, a renewable, GMO-free plant material that can reduce C02 emissions. In fact, their production saves 135 tons in carbon emissions, the equivalent seized from 160 acres of forests in one year.

#### Communicating sustainability

The use of mats made of renewable materials provides additional support in communicating your commitment to environmental sustainability.

#### Recycle

Pharmaclean® mats can be recycled through standard plastic processes.



#### **APPLICATIONS**

- · Pharmaceutical manufacturing
- Microelectronics
- Semiconductors
- Biotech
- Healthcare







#### **Durable and effective**

Adhesive-coated mats offer excellent contamination control performance and resistance to pedestrian and trolley passage.

#### **Competitive** price

Optimal performance combined with material sustainability does not affect the highly competitive cost.







Top protective layer
One-side coated layer
Double-side coated layer
Bottom support layer

TECHNICAL FEATURES				
film material	low density polyethylene			
film thickness	1.4 mil / 35 μm			
adhesive material	acrylic			
adhesive thickness	0.3 mil / 8 μm			
peel adhesion "top"	190 g/25 mm ± 40			
heat resistance	70° C - 48 hours			
efficacy percentage of S. aureus microorganism removal (MRSA)	85%			
efficacy percentage of A. niger microorganism removal	97%			






PRODUCT				
CODE	DESCRIPTION	PACKAGING		
PH707	Pharmaclean® tear-off mat - colors: blue (PH10097), white (PH10099) and grey (PH10098) dimensions cm 61 x 91	8 mats per 30 layers each		
PH706	Pharmaclean® tear-off mat - colors: blu (PH10094), white (PH10096) and grey (PH10095) dimensions cm 46 x 114	8 mats per 30 layers each		
PH708	Pharmaclean® tear-off mat - colors: blu (PH10100), white (PH10102) and grey (PH10101) dimensions cm 66 x 114	8 mats per 30 layers each		
PH709	Pharmaclean® tear-off mat - colors: blu (PH10103), white (PH10105) and grey (PH10104) dimensions cm 91 x 114	8 mats per 30 layers each		
PH710	Pharmaclean® tear-off mat - colors: blu (PH10136), white (PH10137) and grey (PH10138) dimensions cm 46 x 91	8 mats per 30 layers each		
PH713	Pharmaclean® tear-off mat - colors: blu (PH10190), white (PH10191) and grey (PH10192) dimensions cm 46 x 114	4 mats per 60 layers each		
PH715	Pharmaclean® tear-off mat - colors: blu (PH10403), white (PH10404) and grey (PH10405) dimensions cm 91 x 114	4 mats per 60 layers each		



COMPANY WITH QUALITY SYSTEM CERTIFIED BY DNV GL = ISO 9001:2015 =



## STICKY MATS ACCESSORIES



# PHARMACLEAN® CLEANROOM MAT FRAME

#### **DESCRIPTION**

Frame for Pharmaclean® sticky mats, constructed of a hard, rigid, polystyrene base material combined with a non-skid backing. Protect floors from adhesive residue and reduces the risk of slips and falls. Standard white color.



TECHNICAL FEATURES				
material	Polystyrene			
thickness	3,2 mm			
Backing material	Non-Slip polymer coated polyester			

PRODUCTS					
CODE	DESCRIPTION	FOR MATS	PACKAGING		
PH10465	Frame for Pharmaclean® sticky mat, Dimensions 508 x 1194 mm	460 x 1140 mm (PH706 - PH713)	Pack of 1 piece.		
PH10466	Frame for Pharmaclean® sticky mat, Dimensions 660 x 965 mm	610 x 910 mm (PH707)	Pack of 1 piece.		
PH10467	Frame for Pharmaclean® sticky mat, Dimensions 711 x 1194 mm	660 x 1140 mm (PH708)	Pack of 1 piece.		
PH10468	Frame for Pharmaclean® sticky mat, Dimensions 965 x 1194 mm	910 x 1140 mm (PH709)	Pack of 1 piece.		



COMPANY WITH QUALITY SYSTEM CERTIFIED BY DNV GL = ISO 9001:2015 =